

Best practice is to use new N95s. Decontamination does not solve the PPE shortage crisis, and is an emergency practice to be considered during the COVID-19 pandemic. Efficacy and safety of N95 decontamination has not been fully characterized.

COVID N95 DECON & REUSE



UV-C

Use appropriate UV-C source
Validate 1.0 J/cm² dose with sensor
Expose both sides of N95 mask

CORONAVIRUS INACTIVATION

Peer-reviewed data not available for SARS-CoV-2

- +** $\geq 1.0 \text{ J/cm}^2$ of UV-C inactivates* viruses similar to SARS-CoV-2 on N95 FFRs^{1,2**,3}
- $\geq 1.0 \text{ J/cm}^2$ of UV-C yields 2-log reduction of viable *B. subtilis* spores on N95 FFRs⁴
- UV-C light may not reach inner N95 layers for all N95 models⁵
- Elastic straps require additional chemical disinfection¹
- Shadows can block UV-C rays & can leave parts of N95 contaminated

* ≥ 3 -log inactivation

KEY CONSIDERATIONS

- Ensure accurate UV-C dose on all surfaces of N95
- Measure dose at N95 surface with UV-C specific sensor
- Return N95s to original users and ensure handling minimizes cross-contamination
- Perform user seal check before each reuse
- Be aware that data from tests on specific N95 models may not apply to other models

N95 INTEGRITY

- +** N95 keeps fit and filter performance after 10-20 cycles of 1.0–1.2 J/cm² UV-C^{2**,6}
- Each don/doff can reduce N95 fit; some models fit unacceptably after 5 don/doff cycles⁹
- Some damage to N95 seen at high UV-C doses ($\geq 120 \text{ J/cm}^2$)⁶
- Strap and facepiece damage seen on some N95 models after UV-C^{7**,8}

RISKS

- UV light is harmful to eyes and skin;** proper training, engineering controls, and PPE are required before use
- If UV-C source is underpowered, decontamination times may be infeasible
- UV-C may not decontaminate N95 straps or eliminate risk of bacterial co-infection
- Cosmetics and sunscreen on N95 may reduce decontamination efficacy
- Non-uniform irradiance can affect dose, and subsequently, decontamination efficacy

IMPLEMENTATION

- +** Reference documents from University of Nebraska Medical Center⁸ for implementation
- ?** Validate each UV-C source and protocol with a UV-C sensor to ensure adequate dose for decontamination at the N95 surface

CONCLUSION

If implemented properly using sensors to ensure $\geq 1.0 \text{ J/cm}^2$ UV-C dose to the N95, this method likely inactivates SARS-CoV-2; however, this has not yet been confirmed directly with SARS-CoV-2. This method may protect against some bacterial co-infection risks but not all.

SUPPORTING RESEARCH

[1] Mills et al., 2018; [2] Heimbuch & Harnish, 2019**; [3] Lore et al., 2012; [4] Lin et al., 2018; [5] Fisher and Shaffer, 2010; [6] Lindsley et al., 2015; [7] Personal Safety Division, 3M, 2020**; [8] Lowe et al., 2020; [9] Bergman et al., 2012
** = not peer-reviewed

The Content provided by N95DECON is for INFORMATIONAL PURPOSES ONLY and DOES NOT CONSTITUTE THE PROVIDING OF MEDICAL ADVICE and IS NOT INTENDED TO BE A SUBSTITUTE FOR INDEPENDENT PROFESSIONAL MEDICAL JUDGMENT, ADVICE, DIAGNOSIS, OR TREATMENT. Use or reliance on any Content provided by N95DECON is SOLELY AT YOUR OWN RISK. A link to the full N95DECON disclaimer can be found at <https://www.n95decon.org/disclaimer>.



Best practice is to use new N95s. Decontamination does not solve the PPE shortage crisis, and is an emergency practice to be considered during the COVID-19 pandemic. Efficacy and safety of N95 decontamination has not been fully characterized.


COVID N95 DECON & REUSE



UV-C RELATED CONCERNS

UNSUITABLE METHODS

Only **UV-C light with a peak wavelength of 254 nm** has demonstrated substantial germicidal effects on N95 FFRs.^{1**} UV-A (320-400 nm) is not germicidal. UV-B (280-320 nm) has lower germicidal efficiency and has not been validated for N95 FFR decontamination.² **Only use UV-C light sources with a peak wavelength of 254 nm.**



- ✗ Sunlight | Sunlight reaching the Earth's surface does not contain UV-C light;³ there is no evidence in the peer-reviewed literature to support sunlight-assisted decontamination of N95 FFRs.
- ✗ Consumer UV Products | Many consumer UV products do not emit UV-C with sufficient irradiance, and have **peak emission in the UV-A range** (e.g., nail polish curing lamps,⁴ tanning bed lamps,⁵ etc.), which is ineffective for decontamination. Other consumer products may additionally have uniformity or shadowing concerns.
- ✗ < 200 nm UV Sources | UV sources emitting < 240 nm light can produce ozone, which is hazardous to human health.^{6**} Sufficient ventilation is necessary to reduce ozone concentration.⁷
- ✗ Measuring dose from lamp power | UV-C irradiance should not be calculated from rated lamp power, as bulbs do not have 100% efficiency in converting electrical energy to optical power.⁸ Use a **UV-C specific sensor** to measure irradiance at the N95 surface.
- ✗ Doses for air or surface decon | Viral inactivation protocols designed for surfaces or air are insufficient/not effective for N95 decontamination.⁹ Use a substantially higher UV-C dose of 1.0 J/cm² at the N95 surface.^{1**}
- ? Biosafety Cabinets | Many UV-C sources used in research laboratories (e.g., biosafety cabinets) have unacceptable non-uniformity and low power;^{10**} thorough characterization of the UV-C dose at the N95 surface is required for sufficient decontamination.

SUPPORTING RESEARCH

[1] Heimbuch & Harnish, 2019**; [2] Kowalski, 2009; [3] Sagripanti and Lytle, 2005; [4] Dowdy & Sayre, 2013; [5] Spencer & Amonette, 1995
 [6] Oxidation Technologies LLC, 2017;** [7] McClurkin et al., 2013; [8] Lawal et al., 2017; [9] Walker & Ko, 2007; [10] Card et al., 2020**
 ** = not peer-reviewed

The Content provided by N95DECON is for INFORMATIONAL PURPOSES ONLY and DOES NOT CONSTITUTE THE PROVIDING OF MEDICAL ADVICE and IS NOT INTENDED TO BE A SUBSTITUTE FOR INDEPENDENT PROFESSIONAL MEDICAL JUDGMENT, ADVICE, DIAGNOSIS, OR TREATMENT. Use or reliance on any Content provided by N95DECON is SOLELY AT YOUR OWN RISK. A link to the full N95DECON disclaimer can be found at <https://www.n95decon.org/disclaimer>.

